

**DARVALL****We Warm Animals****Before You Start****Using the  
HBC /ZDS****Warm Inspired Air Control Unit**

Use with

**Darvall Heated Smooth-Wall Circuits (HBC)**

On Circle Systems for Cats, Dogs and other small animals up to 80Kgs

or

**Darvall Heated Qubes (ZDS)**

A Non-Rebreathing System for Rodents, Birds, Exotics, Reptiles, Pups and Kittens

**For Veterinary Use Only**

(Not intended, evaluated, or approved for Human Medical purposes)

Certified By MET Laboratories to IEC / EN 61010-1

(Electrical equipment for measurement, control, and laboratory use)

# IMPORTANT

## Before you start using the Warm Inspired Air Control Unit

- Read the Getting Started Guide which is located on the USB drive included with the HBC Controller
- The Controller and Software User Manuals are also on the USB drive
- The Windows HBC Software is on the USB drive and can be installed on your Windows PC to run the monitoring software and download software updates. (Note that the controller will run independently without a computer connection)
- Automated updates of PC Software and Controller Firmware will be available when the Controller is connected to a Windows PC with internet connection by USB cable when the HBC software is started (see Software Manual for more details)
- The website for the latest software and other information is <http://www.aasmedical.net>